

Product Datasheet

NQO1 antibody (orb612235)

Description	NQO1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human NQO1
Target	NQO1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	31 kDa
Uniprot ID	P15559
Dilution Range	WB=1:1000-3000, Flow-Cyt=1:50, ICC=1:100-500, IF=1:50-200
Expiration Date	12 months from date of receipt.