


Product Datasheet

ADAM28 antibody (orb612241)

Description	ADAM28 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ADAM28 (521-620/775 aa)
Target	ADAM28
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	87 kDa
Uniprot ID	Q9UKQ2
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
gastric using
AD...