


Product Datasheet

ALKBH5 antibody (orb612254)

Description	ALKBH5 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human ALKBH5 (1-100/394 aa)
Target	ALKBH5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q6P6C2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Testis (Rat); Lung (Mouse); P...