

Product Datasheet

Parvalbumin antibody [Out of stock] (orb612698)

Description

Mouse monoclonal antibody against

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-594
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Synthetic peptide from the full length Rat Parvalbumin protein
Target	parvalbumin
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:500 dilution of SMC-563 was sufficient for detection of Parvalbumin in 15 ug of Rat Brain by ECL immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2a
Clonality	Monoclonal
MW	12 kDa
Uniprot ID	P02625
NCBI	NP_071944.1
Entrez	25269
Dilution Range	WB (1:500); IHC (1:100)
Expiration Date	12 months from date of receipt.