

Product Datasheet

Alpha Synuclein antibody (pSer129) (orb612713)

Description

Rabbit monoclonal antibody to Alpha

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Human alpha synuclein AA 124-134: AYEMP-pS-EEGYQ-Cys
Target	Alpha Synuclein pSer129
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:500 dilution of SMC-600 was sufficient for detection of Alpha Synuclein pSer129 in 10 µg of Mouse Brain by ECL immunoblot analysis using Goat Anti-Rabbit IgG:HRP as the secondary antibody.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P37840
NCBI	NP_000336.1
Entrez	6622
Dilution Range	WB (1:500)
Expiration Date	12 months from date of receipt.