

Product Datasheet

Alpha Synuclein antibody (pSer129) (orb612717)

Description

Rabbit monoclonal antibody against

conjugated to ATTO-594

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

ATTO-594

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

Human alpha synuclein AA 124-134:
AYEMP-pS-EEGYQ-Cys

Target

Alpha Synuclein pSer129

Preservatives

PBS pH 7.4, 50% glycerol, 0.09%
Sodium azide *Storage buffer may
change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored
at 4°C

Note

For research use only

Application notes

A 1:500 dilution of SMC-600 was
sufficient for detection of Alpha
Synuclein pSer129 in 10 µg of Mouse
Brain by ECL immunoblot analysis using
Goat Anti-Rabbit IgG:HRP as the
secondary antibody.

Isotype

IgG

Clonality

Monoclonal

Uniprot ID
P37840
NCBI
NP_000336.1
Entrez
6622
Dilution Range

WB (1:500)

Expiration Date

12 months from date of receipt.