

Product Datasheet

Tau antibody (pSer202/ pThr205) (orb612734)

Description

Rabbit monoclonal antibody against Tau

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	ATTO-488
Tested Applications	DOT, ELISA, ICC, IHC, WB
Immunogen	Synthetic peptide of Human Phospho Tau (Ser202/Thr205)
Target	Tau (pSer202/ pThr205)
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:500 dilution of SMC-601 was sufficient for detection of Tau in 10 µg Baculovirus by dot blot analysis using Goat Anti-Rabbit IgG:HRP as the secondary antibody.
Isotype	IgG
Clonality	Monoclonal
MW	79 kDa
Uniprot ID	P10636
Entrez	4137
Dilution Range	WB (1:500), ICC/IF (1:500); IHC (1:500)
Expiration Date	12 months from date of receipt.