

Product Datasheet

VPS35 antibody (orb612827)

Description	Mouse monoclonal antibody to VPS35
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Full length recombinant human VSP35
Target	VPS35
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-606 was sufficient for detection of VPS35 in 10 µg of SH-SY5Y by ECL immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary antibody.
Isotype	IgG2b
Clonality	Monoclonal
MW	92 kDa
Uniprot ID	Q96QK1
NCBI	NP_060676.2
Entrez	55737
Dilution Range	WB (1:1000); IP (1:200); IHC (1:100)
Expiration Date	12 months from date of receipt.