

Product Datasheet

Tau antibody (orb612869)

Description	Mouse monoclonal antibody against Tau conjugated to ATTO-594
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	ATTO-594
Tested Applications	DOT, ELISA, ICC, WB
Immunogen	Human Recombinant Tau441 (2N4R), P301S mutant Protein Pre-formed Fibrils
Target	Tau
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1mg/mL
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-608 was sufficient for detection of Tau 2N4R P301S Fibril in 20 ug of Mouse Brain by ECL immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P10636
NCBI	NP_005901.2
Entrez	4137
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.