

Product Datasheet

CRISPR Cas9 Antibody (orb613665)

Description	Our Anti-CRISPR Cas9 primary antibody from Biorbyt is mouse monoclonal. It detects CRISPR Cas9 and i...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant protein from C-terminal region of Cas9 from Staphylococcus aureus.
Target	CRISPR Cas9
Preservatives	PBS + 50% glycerol and 5 mm NaN ₃
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53
Dilution Range	WB: 1:1000, IHC: 1:1000-1:5000, ICC: 1:500-1:5000
Expiration Date	12 months from date of receipt.