

Product Datasheet

GLI1 antibody (orb613759)

Description	Goat polyclonal antibody to GLI1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	GLI1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of U2OS cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of HepG2 cells. Recommended concentration: 10ug/ml.
Protein Sequence	TYSQCPRLEHYGQ
Clonality	Polyclonal
MW	118; 105; 114
Entrez	2735
Dilution Range	P-ELISA: 1:128000,IF: 10 µg/ml, FC: 10ug/ml
Expiration Date	12 months from date of receipt.