

Product Datasheet

CYBB antibody (orb613761)

Description	Goat polyclonal antibody to CYBB
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GP91-PHOX / NOX2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in lysates of cell line Daudi and in Mouse Spleen lysates (calculated MW of 65.3kDa according to Human NP_000388.2 and Mouse NP_031833.3). This band was successfully blocked by incubation with the mmunising peptide, and correspo
Protein Sequence	CDKQEFQDAWKLPK
Clonality	Polyclonal
MW	65.3
Entrez	1536
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.