

## Product Datasheet

### GITR antibody (orb613834)

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit monoclonal antibody to GITR                                                      |
| <b>Species/Host</b>        | Rat                                                                                     |
| <b>Reactivity</b>          | Mouse                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                            |
| <b>Tested Applications</b> | FC                                                                                      |
| <b>Immunogen</b>           | Mouse GITR.                                                                             |
| <b>Target</b>              | GITR                                                                                    |
| <b>Preservatives</b>       | PBS with 0.02% Proclin 300.                                                             |
| <b>Concentration</b>       | 1 mg/ml                                                                                 |
| <b>Storage</b>             | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.        |
| <b>Note</b>                | For research use only                                                                   |
| <b>Application notes</b>   | This antibody binds to murine GITR and can be used to stain for GITR in flow-cytometry. |
| <b>Isotype</b>             | IgG                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                              |
| <b>Purity</b>              | Purified                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                         |