

Product Datasheet

A2A-Adenosine R antibody (orb613835)

Description	Mouse monoclonal antibody to A2A-Adenosine R
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA
Immunogen	Human recombinant A2-adenosine receptor.
Target	A2A-Adenosine Receptor
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to the A2 adenosine receptor, a G-protein-coupled receptor for adenosine that, upon ligand binding, activates adenylyl cyclase. The antibody acts as an allosteric inhibitor.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P29274
Expiration Date	12 months from date of receipt.