

Product Datasheet

GITR antibody (orb613890)

Description	Rabbit monoclonal antibody to GITR
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	nB2/6TG cells transfected with mouse GITR injected into DA rats.
Target	GITR
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	A chimeric version of this antibody featuring a human IgG1 Fc-portion engineered to lack the N-linked glycosylation site at amino acid 297 was shown to enhance humoral immune responses by blocking GITR signalling (Ponte 2010).
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.