

Product Datasheet

Glycophorin A M antigen antibody (orb613929)

Description	Rabbit monoclonal antibody to Glycophorin A M antigen
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	M2A1 was prepared by immunizing BALB/c mice with human red blood cells of the M phenotype, and was screened by testing for anti-M activity.
Target	Glycophorin A M antigen
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	M2A1 can bind to M antigen of glycophorin A and can be used in ELISAs, immunoblots and modified M glycoprotein can be used for inhibition tests.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P02724
Expiration Date	12 months from date of receipt.