



## Product Datasheet

### CD74 antibody (orb613943)

<b>Description</b>	Rat monoclonal antibody to CD74
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IHC-Fr, IP
<b>Immunogen</b>	Mouse lymphocytes of strain A injected into Wistar rats.
<b>Target</b>	CD74
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	The In-1 antibody was used to identify CD74 on lymphocytes in mice. It can be used in immunoprecipitation as well as in flow cytometry (Takahashi et al. 2009, PMID: 19413900) and frozen-section IHC (Momburg et al. 1986, PMID: 3079804).
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	CD74
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Mouse spleen (A...)



Immunofluorescence analysis of paraforma...