

Product Datasheet

DARC antibody (orb613953)

Description	Rabbit monoclonal antibody to DARC
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC
Immunogen	CHO cells transfected with cDNA encoding the Duffy glycoprotein.
Target	DARC
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The antibody binds to DARC, a glycoprotein chemokine receptor and the only known receptor for Plasmodium vivax.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q4VBN9
Expiration Date	12 months from date of receipt.