


Product Datasheet

Bovine IGF-1 ELISA Kit (orb614011)

Description	Bovine IGF-1 ELISA Kit
Reactivity	Bovine
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G49-A118
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	21841 MW
Uniprot ID	P05019
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Bovine IGF-1
PicoKine
ELISA Kit
standard...