

## Product Datasheet

### CEACAM5 Antibody (orb614107)

## Description

CEACAM5 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested Applications

ELISA, IHC, WB

### Immunogen

E.coli-derived human CEACAM5 recombinant protein (Position: K35-N650).

### Form/Appearance

Lyophilized

### Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

### Note

For research use only

### Application notes

"Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users." . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

### Isotype

Rabbit IgG

### Clonality

Polyclonal

### MW

130-200 kDa

### Uniprot ID


**P06731**

### Dilution Range


"Western blot, 0.25-0.5µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat Direct ELISA, 0.1-0.5µg/ml, Human"

### Expiration Date


12 months from date of receipt.



IHC analysis of CEACAM5 using anti-CEACA...



IHC analysis of CEACAM5 using anti-CEACA...



IHC analysis of CEACAM5 using anti-CEACA...