

Product Datasheet

DDX1 Antibody (orb614109)

Description

DDX1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications

Immunogen

E.coli-derived human DDX1 recombinant protein (Position: K562-F740).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

"Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users." . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

88 kDa

Uniprot ID


[Q92499](#)

Dilution Range


"Western blot, 0.1-0.25µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human
Flow Cytometry, 1-3µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5µg/ml, Human"

Expiration Date


12 months from date of receipt.



Flow
Cytometry
analysis of
SiHa cells us...



IF analysis of
DDX1 using
anti-DDX1
anti...



IHC analysis of
DDX1 using
anti-DDX1
ant...