

Product Datasheet

DHX15/prp43 Antibody (orb614111)

Description

DHX15/prp43 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

A synthetic peptide corresponding to a sequence of human DHX15/prp43 (DIKPEWLVKIAPQYYDMSNFPQCEAKRQLDRIIAKLQSKEYSQY).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

91 kDa

Uniprot ID


O43143

Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml,
Human Immunocytochemistry/Immunofluorescence, 5µg/ml,
Human Flow Cytometry(Fixed), 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



WB analysis using anti-DHX15/prp43 antibody...



IF analysis of DHX15/prp43 using anti-DH...



Flow Cytometry analysis of Ana-1 cells u...