

Product Datasheet

TRAF6 Antibody (orb614140)

Description
TRAF6 Antibody
Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IHC, WB

Applications
Immunogen

E.coli-derived mouse TRAF6 recombinant protein (Position: Q14-D512).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

"Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users." . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


60-65 kDa

Uniprot ID
P70196
Dilution Range


"Western blot, 0.25-0.5µg/ml, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat Flow Cytometry, 1-3µg/1x10⁶ cells, Mouse Direct ELISA, 0.1-0.5µg/ml, Mouse"

Expiration Date


12 months from date of receipt.



WB analysis using anti-TRAF6 antibody.La...



IHC analysis of TRAF6 using anti-TRAF6 a...



IHC analysis of TRAF6 using anti-TRAF6 a...