

Product Datasheet

HER2 (phospho-Tyr1248) antibody (orb6146)

Description	HER2 (phospho-Tyr1248) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human HER2 around the phosphorylation site of Tyr1248 (PE(p-Y)LG)
Target	Phospho-HER2 (Tyr1248)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	138 kDa
Uniprot ID	P04626
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.