

Product Datasheet

CaMKII (Phospho-Thr 286) Antibody (orb615040)

Description	CaMKII (Phospho-Thr 286) Antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q9UQM7
Expiration Date	12 months from date of receipt.