

Product Datasheet

HER2 / ErbB2 (Phospho-Thr 686) Antibody (orb615046)

Description	HER2 / ErbB2 (Phospho-Thr 686) Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P04626
Expiration Date	12 months from date of receipt.