

Product Datasheet

c-Jun(Phospho-S63) Antibody (orb615079)

Description	c-Jun(Phospho-S63) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05412
Expiration Date	12 months from date of receipt.