

Product Datasheet

Acetyl and phospho-Histone H3 (Ac-K9/p-Ser10) antibody (orb6162)

Description Acetyl and phospho-Histone H3 (Ac-K9/p-Ser10) antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen KLH conjugated synthetic peptide derived from human Histone H3 around the site of Acetyl-Lys9 and phospho-Ser10 (AR(Ac-K)(p-S)TG)

Target Histone H3 (acetyl K9, phospho S10)

Preservatives 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance Liquid

Concentration 1mg/ml

Storage Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal


MW 15 kDa

CAS Number Histone H3.1

Uniprot ID **P68431**

Dilution Range WB=1:500-2000

Expiration Date 12 months from date of receipt.



Western blot: Lane 1: Pancreas (Rat); Lane ...