

Product Datasheet

[KO Validated] MAP2K2 Polyclonal Antibody (orb616352)

Description	[KO Validated] MAP2K2 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P36507
Expiration Date	12 months from date of receipt.