

Product Datasheet

Acetyl-Histone H2A (K5) antibody (orb6164)

Description	Acetyl-Histone H2A (K5) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human Histone H2A around the acetylation site of Lys5 (RG(Ac-K)QG)
Target	Histone H2A (Acetyl K5)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14 kDa
Uniprot ID	P0C0S8
Expiration Date	12 months from date of receipt.