

## Product Datasheet

### [KO Validated] BUB1 Polyclonal Antibody [Out of stock] (orb617469)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | [KO Validated] BUB1<br>Polyclonal Antibody |
| <b>Species/Host</b>        | Rabbit                                     |
| <b>Tested Applications</b> | IHC-P, WB                                  |
| <b>Storage</b>             | Store at -20°C                             |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Polyclonal                                 |
| <b>Uniprot ID</b>          | <b>O43683</b>                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.            |