



Product Datasheet

cTnT antibody (orb6180)

Description	cTnT antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Troponin T, cardiac muscle. (201-298/298 aa)
Target	Cardiac Troponin T
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	P45379
Dilution Range	Flow-Cyt=1µg /Test, IHC-P=1:100-500, IHC-F=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of Thymus Cell ...



Immunohistochemical staining of Rat hear...