



Product Datasheet

IL26 antibody (orb6198)

Description	IL26 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human IL26. Please contact us for the exact immunogen sequence. The peptide is available as orb12919.
Target	IL26
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	20 kDa
Uniprot ID	Q9NPH9
Entrez	55801
Dilution Range	IHC-P: 1:100-500, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical staining of rat thyr...



Western blot analysis of rat heart (lane...