

Product Datasheet

Alpha-tubulin Antibody (orb620993)

Description	Alpha-tubulin Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P68366
Expiration Date	12 months from date of receipt.