



Product Datasheet

IL6 antibody (orb6210)


Description	IL6 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human IL6. Please contact us for the exact immunogen sequence. The peptide is available as orb374783.
Target	IL6
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	24 kDa
Uniprot ID	P20607 , P08505
NCBI	031168 , 112445 , 11
Entrez	16193
Dilution Range	IF/ICC: 1:100-800, WB:1:200-1000, IHC-P:1:100-800
Expiration Date	12 months from date of receipt.



IF image of mouse skin tissue using anti...



IHC-P staining of mouse skin tissue usin...



Western blot analysis of mouse liver (la...