


Product Datasheet

SNAP-25 Antibody (orb621147)

Description	SNAP-25 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P60880
Dilution Range	WB: 1:1000, IF/ICC: 1:100, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis on SH-SY-5Y cells ...



Immunofluorescent staining of SNAP25 in ...