



Product Datasheet

CD11c / ITGAX Antibody (orb621157)

Description	CD11c / ITGAX Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P20702
Expiration Date	12 months from date of receipt.



Western blot analysis of Integrin αX ex...



Western blot analysis of Integrin αX ex...