



Product Datasheet

IL7 Receptor alpha (phospho-Tyr449) antibody (orb6212)

Description

IL7 Receptor alpha (phospho-Tyr449)


Immunohistochemical
analysis of paraffin...

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IL-7Ra around the phosphorylation site of Tyr449 (EA(p-Y)VT)
Target	phospho-IL7R (Tyr449)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	P16871
Dilution Range	WB=1:500-2000, Flow-Cyt=2ug/Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.