

Product Datasheet

Integrin beta 1 Antibody (orb621464)

Description	Integrin beta 1 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05556
Expiration Date	12 months from date of receipt.