




## Product Datasheet

### Myoglobin Antibody (orb621770)

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | Myoglobin Antibody                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                    |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | ICC, IHC, IP, WB                                                                          |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                     |
| <b>Clonality</b>           | Recombinant                                                                               |
| <b>Uniprot ID</b>          | <b>P02144</b>                                                                             |
| <b>Dilution Range</b>      | WB: 1:1,000-5,000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western blot analysis of Myoglobin on hu...



Immunohistochemical analysis of paraffin...