

Product Datasheet

Insulin Receptor Beta (phospho-Tyr1345) antibody (orb6219)

Description

Insulin Receptor Beta (phospho-Tyr1345)


Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunohistochemical
analysis of formalin...

Immunogen

KLH conjugated Synthesised
phosphopeptide derived from human
Insulin Receptor alpha around the
phosphorylation site of Tyr1355 (RS(p-
Y)EE)

Target

Phospho-Insulin Receptor Beta (Tyr1355)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,
0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one
year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

68/152 kDa

CAS Number

3

Uniprot ID
P06213
Dilution Range

WB=1:500-2000, IHC-P=1:100-500, Flow-
Cyt=3ug/Test, IHC-F=1:100-500, IF=1:100-
500, ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.