



Product Datasheet

Thymidine Kinase 1 Antibody (orb622100)


Description	Thymidine Kinase 1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P04183
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T(Lane 1), H...



Immunohistochemical analysis of paraffin...



Immunofluorescence staining of HeLa cell...