



Product Datasheet

DDIT3 Antibody (orb622172)


Description	DDIT3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P35638
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-12 cells lys...



Immunohistochemical analysis of paraffin...



Flow cytometric analysis of Hela cells
W...