

Product Datasheet

DNA Polymerase beta Antibody (orb622357)

Description	DNA Polymerase beta Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P06746
Expiration Date	12 months from date of receipt.