

Product Datasheet

Estrogen Inducible Protein pS2 Antibody (orb622425)

Description	Estrogen Inducible Protein pS2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P04155
Expiration Date	12 months from date of receipt.