

Product Datasheet

SMAD2 (Phospho-Ser465/pSer467)/SMAD3 (Phospho-Ser423/pSer425) antibody [Out of stock] (orb622944)

Description	Rabbit polyclonal antibody to SMAD2 (Phospho-Ser465/pSer467)/SMAD3 (Phospho-Ser423/pSer425)
Species/Host	Rabbit
Tested Applications	WB
Target	SMAD2 (pSer465/pSer467)/SMAD3 (pSer423/pSer425)
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.