

Product Datasheet

IL-10 antibody (orb623012)

Description	Rat monoclonal antibody to IL-10
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-P, WB
Immunogen	Human IL-10 expressed by COS cells.
Target	IL-10
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.