

Product Datasheet

CD8 ALPHA antibody (RPE) (orb623027)

Description	Mouse monoclonal antibody against CD8 ALPHA conjugated to RPE
Species/Host	Mouse
Reactivity	Porcine
Conjugation	RPE
Tested Applications	FC
Immunogen	Porcine lymphocytes
Target	CD8 ALPHA
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.