

Product Datasheet

CD107a antibody (APC) (orb623046)

Description	Mouse monoclonal antibody against CD107a conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Human PBMCs
Target	CD107a
Preservatives	0.02% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Allophycocyanin (APC) - liquid
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.