

Product Datasheet

CD11a antibody (orb623053)

Description	Mouse monoclonal antibody to CD11a
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IP
Immunogen	Leucocytes from a patient suffering from a LGL-type leukemia
Target	CD11a
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.